



HugeMed

## Single-use Ureterorenoscope

**6.3Fr** **3.6Fr**

Outer Diameter

Working Channel

3 o'clock & 9 o'clock



9 o'clock  
exit



3 o'clock  
exit



### Enhanced Bending Section

285° up and down bending angle offers much compensation even with placement of instruments or laser fibers.



### Insertion Tube Rotation

Easier to explore and facilitate team collaboration



### High Image Quality

Provides exceptional visibility within the urinary tract, and delivers consistent, high-quality illumination throughout the procedure.



## Single-use Ureterorenoscope

### Specifications




#### Single-use ureterorenoscope-HU30M

Insertion Tube Diameter	6.3Fr	Working Length	700 mm
Working Channel Diameter	3.6Fr	Total Length	905 mm
Bending Angle	285° up/down	View Angle	120°
Insertion Rotation	120° left/right	View Depth	3-50 mm

# Single-use Ureteral Access Sheath


- Improve surgical field of vision
- Reduce RIRS operating time
- Increase stone-free rate in a single-session

Available in  
Multiple Sizes



## Dilator tip

The dilator features excellent flexibility, with a smooth, tapered tip and seamless transition to the sheath, facilitating dilation.




## Slider Valve

Sliding Valve adjusts vent opening for precise intraoperative negative pressure control.



## Bending Part

Ultra-flexible distal segment with stainless steel coil support, enabling large-angle deflection with the endoscope.



## Snap-fit

Guarantees coaxial advancement with zero separation risk. One-hand quick unlock for efficient, precise access.



## Stainless steel coil-reinforced sheath

Stainless steel coil reinforcement for robust support, high collapse resistance under negative pressure, ensuring safe procedures.

**Shenzhen HugeMed Medical Technical Development Co., Ltd.**  
1601, 1701, 1801, Luohu New Materials Building, No. 6, Qingyuan Road, Qingshuihe Community,  
Qingshuihe Street, Luohu District, Shenzhen, China, 518000



+86-755-22275833



+86-755-22275833



info@hugemed.net



hugemed.net

# HugeMed